

DELUXE ADJUSTABLE FIT COOLING VEST (#KMVZ)

WITH (5 TO 12) 4.5" X 6" COOLING PACKS



THIS VEST IS COMPATIBLE WITH BOTH KOOL MAX® AND COOL58® COOLING PACKS.

STYLE CHOSEN:

ZIPPER VEST PONCHO VEST VELCRO VEST

MATERIAL CHOSEN:

COTTON BREATHABLE, EASY-CLEAN NYLON

COOLING PACKS CHOSEN:

KOOL MAX® WATER-BASED COOLING PACKS

Reusable, water-based packs require access to a freezer and offer the highest level of cooling energy. It's science! It takes more energy to charge these packs in a freezer and in return, the packs offer greater cooling energy for a longer time.

COOL58® PHASE CHANGE COOLING PACKS

Reusable, phase change packs fully charge (freeze) at 58°F (14.4°C) and provide a comfortable, constant cooling temperature of 58°F. Ideal choice when freezers are not available; these packs can be activated in 20 minutes in a container of ice water! Additional time required if activating packs in the refrigerator or freezer.

MINIMAL WEIGHT, MAXIMUM COVERAGE

The number of pack pockets sewn into the vest varies according to vest size. This allows maximum cooling pack coverage according to the size of the person while minimizing the weight. The location of the packs along the spine in the back also maximizes the cooling effect while minimizing weight.

FLEXIBILITY TO MEET YOUR COOLING NEEDS

Individual packs allow maximum flexibility in varying the number of packs used. This allows for customizing the vest weight and location of cooling.

QUALITY MATERIALS

Polar's cooling vests are designed and constructed for long term performance. They will last for many years with proper care, maintenance and storage.

IMPROVES PERFORMANCE, SAFETY AND COMFORT

Cooling lessens the debilitating effects of heat during work, athletics and day-to-day life.



If you have any questions or concerns regarding your product, please contact our customer service at 800.763.8423 or send an e-mail to polar@polarproducts.com. We will be happy to help!

ACTIVATION & USE INSTRUCTIONS

WARNINGS AND CAUTIONS

KOOL MAX® NON-TOXIC WATER-BASED COOLING PACKS:

Polar Products recommends that all cooling products activated in a freezer be worn over clothing or other fabric layers to regulate the products' cooling potential.

The amount of layering required to ensure the highest level of cooling and comfort will vary from person to person and is strictly the responsibility of each individual user.

COOL58® NON-TOXIC COOLING PACKS:

In case of pack rupture or leakage, discontinue use immediately. Although encapsulated material is non-toxic, it can cause a skin rash. Flush exposed area with warm, soapy water and wash with baking soda.

Contact physician immediately if swallowed or in case of eye exposure. Wash all exposed products with baking soda and detergent in warm water.

Do not expose cooling packs to heated equipment and/or environments exceeding 300°F.

Note: Please contact us to request a Material Safety Data Sheet.

CAUTION:

Consult your healthcare provider prior to using this or any cooling products for a medical condition.

Polar Products, Inc. will not be liable for injuries or damages resulting from misuse or misapplication of our products.

WEARING YOUR VEST

Our high-quality vests are constructed of cotton or nylon fabric. The vest should fit snug to the body for optimal performance. The vest provides adjustable fit through the use of hook and loop straps at the shoulders, chest and waist. Simply place the vest on your body and adjust for a snug fit.



PLACEMENT OF PACKS

The vest uses 5 to 12 individual 4.5" x 6" cooling packs. Not all packs must be used with the vest. Each person may customize the number of packs used and the location of the packs. Some people find they may only need a few packs on the back for instance, others want to use them all. It is your choice!

USING THE KOOL MAX® COOLING PACKS

Place the packs flat into a freezer until frozen hard. For best results, store the packs in the freezer, so they're ready when you are. When ready to use, simply remove the packs from the freezer and place in the insulated pockets in the vest. Frozen packs may be transported or stored in a cooler for later use.

USING THE COOL58® COOLING PACKS

Place the packs flat into a container of ice water, refrigerator or freezer until fully-charged (frozen hard). The packs turn solid white when they're ready. When ready to use, simply remove the packs from the ice water, refrigerator or freezer and place in the insulated pockets in the vest.

Note: Additional time required if activating packs in the refrigerator or freezer.

COOLING DURATION

Kool Max® packs provide cooling for up to 4 hours in any climate; Cool58® packs provide cooling for up to 3 hours in any climate. Cooling times are estimates.

The length of time Kool Max® and Cool58® cooling vests will cool depends on many factors, including the activity and body metabolism of the person wearing the vest and the temperature of the environment. The cooling potential may also be adjusted by wearing one or more layers of clothing between the body and the vest.

Note: One or more layers of insulating fabric should always be worn between the skin and the cooling pack when using the packs from the freezer.

CARE AND MAINTENANCE

Vest: For cotton vests, hand wash in cool water with a mild detergent. Hang in a well-ventilated area to dry. For nylon vests, wipe down with a damp cloth or anti-bacterial disinfecting wipe. It may also be cleaned by hand washing in cool water with a mild detergent. Hang in a well-ventilated area to dry.

Cooling Packs: Packs may be wiped clean with a damp cloth.

Note: Store packs in freezer or cool, dry environment when not in use.

TESTIMONIALS

HEAR WHAT OUR CUSTOMERS HAVE TO SAY ABOUT OUR COOLING VESTS!

"Your cooling vests are amazing. I have RA and take Methotrexate, and as a result I am very heat intolerant. I received my vest last week and was able to be outside in 100-degree heat for two hours and didn't over heat! Thank you for making these life changing products."

– Holli in Texas

"My first time wearing it was in the Midwest's hot and humid weather. It kept me so comfortably cool that I forgot I was wearing it. I am so happy to have this vest! I can do things that the heat would make impossible to do without the vest. A life changer!"

– Judy in Minnesota

"Living in Florida, I can't live without this. I tell people about it all the time, and when people ask me what I'm wearing, I make them try it on. They are always so impressed and wanting one for themselves! This is the most amazing product invented on earth."

– Margo in Florida

"My wife has had a Polar Kool Max vest for over 10 years to help her deal with the heat - she has MS. We're looking at five or six days with temps around 90 this week, and she would be stuck in the house without this cooling vest. It is invaluable at times like this."

– Chip in California

"My first experience with the vest yesterday was wonderful. Absolutely met my expectations and the craftsmanship of the vest and subsequent cooling pack pockets is really well done. Thanks for making a great product. I purchased it after much research and I believe I made the right decision."

– Johnathan in North Carolina

ADDITIONAL PRODUCTS

POLAR OFFERS SOLUTIONS FOR COOLING EVERY BODY®!

There is no "one size fits all" for cooling vests, and we have a variety of cooling products allowing you to customize your individual cooling needs.



WOMEN'S FASHION VESTS



MEN'S FASHION VESTS



HATS



TORSO VESTS



NECK WRAPS



HEAD WRAPS



ANKLE & WRIST WRAPS



BREAST & BRA COOLERS

Polar Products provides natural, drug-free relief from the heat! Visit polarproducts.com to see our full line of body cooling technologies and products. For customer service, call 800.763.8423.



SAVE MONEY WITH A COOLING KIT!

WE'VE BUNDLED OUR MOST POPULAR COOLING VESTS AND ACCESSORIES INTO KITS SOLD AT A DISCOUNTED PRICE!

Use the products separately as needed or together for maximum cooling.



HOME & AWAY KITS:

Our "Home and Away" Kits include both cooling pack technologies: a set of Kool Max® frozen water-based packs and a set of Cool58® 58°F phase change packs. Kool Max® packs are the best choice when you're at home or anywhere with access to a freezer. Cool58® packs can be reactivated on-the-go in a cooler of ice water when you're traveling!

These kits are available with every style of cooling pack vest.



COOL58® PHASE CHANGE COOLING ACCESSORIES KIT:

Includes our best cooling accessories for before, during, and after exercise: your choice of Cool58® Hat, a Cool58® Deluxe Neck Wrap, and a pair of Cool58® Performance Quick Dry Wrist Wraps. Accessories target key arteries for the ultimate cooling.



COOL COMFORT® SPORTS KIT:

Includes a Sports Vest, a neck band, and a pair of wrist and ankle wraps. This hybrid cooling system combines the efficient cooling power of evaporation with the direct cooling of frozen water.

ALSO AVAILABLE FROM POLAR:

HOT/COLD THERAPY PRODUCTS



ACTIVE ICE® COLD THERAPY SYSTEMS

The Active Ice® circulating cold water therapy system is great for drug-free pain relief after injury or post-surgery. With more than ten different body-part specific therapy pads, we have you covered!

These systems are designed to be used anywhere and provide effective treatment for reducing pain and inflammation, and facilitating recovery. Great for knees, shoulders, backs, hips, ankles, elbows, calves, limbs, and more.



SOFT ICE®

Reusable cold/hot packs and compression wraps to reduce inflammation.



THERA-TEMP®

Microwavable moist heat wraps offer relief of chronic pain.



Polar Products is a family-owned company with over 35 years of experience manufacturing effective, high-quality and affordable body cooling and therapy products.

polarproducts.com
800.763.8423

