



Fashion Scarves

Reusable Quick Cool™ packs (for short-term relief) or Kool Max® cooling packs (for extended relief) can be used in one comfortable, fashionable scarf for discreet cooling.



How does it work?

There are two pack options:

- **Quick Cool™ cooling packs** stay soft when frozen for comfortable, short-term cooling for up to 45 minutes.
- **Kool Max® water-based packs** freeze solid in the freezer and offer extended cooling for up to 2 hours.
- Both pack technologies fit into pockets sewn in the scarf, which ties comfortably around the head or neck.

CARE & USE MANUAL



The Advantages of Quick Cool™ and Kool Max®:

Discreet Relief

Small, lightweight cooling packs tuck discreetly in pockets sewn into fashionable cotton scarves!

Quality Materials

Fashion scarves are designed and constructed for long-term performance. They will last for many years with proper care, maintenance and storage.

Highest Cooling Energy

Kool Max® and Quick Cool™ water-based technology provides the highest cooling energy of any cooling pack. These packs require the use of a freezer but return this energy in cooling. It's science!

Improves Performance, Safety and Comfort

Cooling lessens the debilitating effect of heat during work, athletics and day-to-day life for those with heat sensitivity.

Natural, Drug-Free Cold Therapy

Not only will Quick Cool™ packs keep you cool... cold therapy is excellent drug-free pain relief for bruises, muscle aches, swelling, headaches, chronic pain and much more! The scarf can be used to position the pack on the target area.

Caution: Use only the Quick Cool™ packs - not the Kool Max® packs - for therapy applications. Always layer insulating fabric between the pack and the skin. Cold therapy should be limited to 20 minute intervals. Consult your health care professional.



Activation and Use Instructions

If you have any questions or concerns regarding your product, please contact our Customer Service at **1.800.763.8423** or send an e-mail to polar@polarproducts.com. We will be happy to help!

Using the Cooling Packs

Our reusable Quick Cool™ and Kool Max® Packs are water-based, non-toxic and non-flammable. Place the packs individually and flat into a freezer. For best results, store the packs in the freezer, so they're ready when you are.

Note: Quick Cool™ Packs will remain soft and flexible when frozen. Kool Max® Packs will freeze hard and solid.

When ready to use, simply remove the packs from the freezer and place in the pockets sewn in the scarf. Frozen packs may be transported or stored in a cooler for later use.

Wearing the Fashion Scarves

Scarves should fit snug to the body for optimal performance. The scarves can be tied around the neck or head for a comfortable, stylish fit.

Temperature Adjustment

The cooling potential may also be adjusted by wearing one or more layers of clothing between the body and the garment. **Note: one or more layers of insulating fabric should always be worn between the skin and the cooling packs.**

Maintenance and Care

Fashion Scarves: Hand wash in cool water with a mild detergent. Hang in a well-ventilated area to dry.

Cooling Packs: Packs may be wiped clean with a damp cloth. **Note: Store packs in freezer or cool, dry environment when not in use.**

Cooling Duration

The length of time fashion scarves will cool depends on many factors, including the activity and body temperature of the person wearing the scarf and the temperature of the environment. **Quick Cool™ Packs are soft and flexible when frozen and will therefore thaw more quickly than Kool Max® Packs, which are frozen solid. Quick Cool® Packs are intended for short-term relief, typically up to 45 minutes. Use the Kool Max® Packs for extended cooling relief for up to two hours.**

Quick Cool™ Packs



Kool Max® Packs



Cautions and Warnings: Consult your health care professional prior to using this or any cooling products for a medical condition. Kool Max® and Quick Cool™ cooling packs are designed for the specific objective of cooling the body. Polar Products recommends that all cooling products be worn over clothing or other fabric layers to regulate the products' cooling potential. The amount of layering required to ensure the highest level of cooling and comfort will vary from person to person and is strictly the responsibility of each individual user. **Polar Products Inc. will not be liable for injuries or damages resulting from misuse or misapplication of our products.**