

Active Ice® System

- Convenient and effective cold therapy and compression system
- Your choice of therapy bladder for optimal targeted therapy
- Optional compression wraps hold the universal bladder comfortably on the targeted body part and insulate the area
- Y-shaped hose connector allows for multi-wrap or multi-person simultaneous use

Active Ice® Head Cap

Provides the highest level of cold or heat therapy for relief of migraines. Two styles of full coverage caps offer continuous and comfortable therapy.



Active Ice® Head Cap

POLAR Products Active Ice® System



Also available from POLAR Products:

Soft Ice® Cold/Hot Packs and Compression Wraps



SOFT AND FLEXIBLE WHEN FROZEN

Thera-Temp® Moist Heat Therapy Wraps



DEEP-PENETRATING MOIST HEAT

"I call it my portable life support system - it provides immediate, continuous pain relief that is non-sedating and not addictive. I love this machine and it goes everywhere with me - day and night."

- Malinda, a physician in Manistee, MI

#AIS-D



Active Ice® Double Therapy System

offers convenient cold therapy to two parts of the body simultaneously

Custom systems available!

Contact us at 1.800.763.8423 or polar@polarproducts.com for a quote.

Visit www.polarproducts.com to see our complete hot and cold therapy product line.



GSA Contract
#V797D-30220

Therapy Bladders

Universal U-Shaped Bladder can be used on nearly any part of the body with the included 3" x 44" elastic belt.



#UB
U-Shaped Bladder



#RB
Rectangular Bladder



#KAW
Knee



#LAW
Lumbar



#FB
Foot Bladder

Targeted Compression Wraps

Optional wraps hold the U-shaped bladder comfortably to the body, provide a gentle compression and insulate the targeted area.



#SAW
Shoulder

Active Ice® System

Circulating Cold Water Therapy System for Post-Surgery Recovery

9-Quart or 15-Quart Cooling Reservoir offers convenient extended cooling and minimizes ice refill. Easy to use for effective cold therapy; simply fill the reservoir with ice and water!

9-quart cooling reservoir shown on front cover

"I truly think this is one of the best pieces of equipment for reducing swelling that I have ever used." - Greg in Phoenix, AZ

Adjustable Flow Rate Temperature Control

30-Minute Timer, On/Off Switch or Both

Power Options:

- Lithium Ion Battery Pack
- 12 volt DC Adapter

System includes plug-in grounded transformer



Retractable Handle and Wheels on 15-quart cooler allow for mobility

High-Flow, Dry, Quick Disconnect Couplings allow freedom of movement if needed

5-Feet or 8-Feet of Insulated Water Lines circulate cold water to the therapy bladder



www.polarproducts.com • 1.800.763.8423

