

# POLAR PRODUCTS SKULL CAP ACCESSORY FOR POLAR COOLING HEAD CAP



## CAN BE USED IN TWO WAYS:

1. As a barrier between the cooling head cap and bare skin when there is no hair or very short hair - to give protection to the skin.
2. Alternatively, the skull cap may be used over top of the Polar cooling head cap as a way to get better contact between the cooling cap and the head for better cooling of the scalp.



**NEED ASSISTANCE? CALL US - WE'RE HAPPY TO HELP!**  
**polarproducts.com | 800.763.8423 | polar@polarproducts.com**

*If you need an extra skull cap to protect the skin or wear over the cooling cap, call us—we'll send one free of charge.*