

# CHOOSING THE IDEAL COOLING PACK BASED VEST OR ACCESSORY



## 1 WHICH REUSABLE COOLING TECHNOLOGY SHOULD YOU CHOOSE?



### KOOL MAX®

Cooling vests with Kool Max® frozen water-based cooling packs fit into insulated pockets. Packs require access to a freezer and 2 to 3 hours to charge. They offer the highest level of cooling energy. *It's science! It takes more energy to charge these packs in a freezer and in return, the packs offer greater cooling energy for a longer time.* Vests will cool for approximately 3 to 4 hours; accessories will cool for 1 to 2 hours.



### COOL58®

Cooling vests with Cool 58® phase change cooling packs fit into insulated pockets. Packs freeze at 58°F (14.4°C) and provide a comfortable, constant cooling temperature of 58°F. Ideal choice when freezers are not available; these packs can be recharged in a cooler of ice water in just 15 to 20 minutes or a refrigerator in 1-1.5 hours. Vests will cool for approximately 2 to 3 hours; accessories will cool for 1 to 1.5 hours.



### HOME AND AWAY KITS

Our home and away kits include a vest and one set of each style of pack. Kool Max® for when you have easy access to a freezer; Cool 58® when you are traveling and only have access to a bucket of ice water.

## 2 WHICH STYLE OF VEST WOULD YOU PREFER?




### FITTED FASHION VESTS

This vest comes in both a women's and men's style. It has a clean look with no adjustments and is lightweight with 5 to 10 cooling packs depending on size. It can be worn over or under clothing. The fashion vest comes in a 97% cotton / 3% spandex blend that provides a snug fit with a little stretch.



### ADJUSTABLE VESTS

This cooling vest has adjustable hook and loop straps at the chest, waist and shoulders for a snug, comfortable and customized fit.





### WOMEN'S FASHION VESTS

This fitted vest includes 5 to 8 packs depending on size. This configuration also allows the user to adjust the number, location and weight of the vest as desired. Sizes range from XS to 4XL. Each pack weighs 6.4 oz. Back cooling packs are located along the spine for efficient cooling.

The women's vest comes in black, khaki, light blue, and pink. There are two outside side pockets for general use.



### MEN'S FASHION VESTS

This fitted vest includes 8 to 10 packs depending on size. This configuration also allows the user to adjust the number, location and weight of the vest as desired. Sizes range from S to 4XL. Each pack weighs 6.4 oz. Back cooling packs are located along the spine for efficient cooling.

The men's vest comes in black, gray and khaki. There are two outside side pockets for general use.



### ADJUSTABLE VESTS

Adjustable vests include 8 to 12 packs depending on vest size. This configuration allows the user to adjust the number, location and weight of the vest as desired. Sizes range from XS/S to 3XL/4XL. Each pack weighs 6.4 oz. Back cooling packs are located along the spine for efficient cooling.

Adjustable vests are available in khaki, blue, gray and black.

**"Polar cooling vests are amazing. I have RA and take Methotrexate, and as a result I am very heat intolerant. I received my vest last week and am now able to be outside in 100-degree heat for two hours and didn't over heat! Thank you for making these life changing products."**



**"My daughter has Mitochondrial Disease, and she is unable to regulate her body temperature. This vest has given her the freedom to go outside when it is over 75 degrees. THANK YOU for this precious gift, and helping her have the freedom to go places and still maintain a safe core temp."**



**"Polar Fashion Cooling Scarves have really helped me deal with MS-related heat intolerance. I have them in Black, Khaki and Blue. They are fashionable enough that people don't think that they are a medical aid! I have worn them in the sweltering August heat while traveling throughout Asia and Europe. They have been a tremendous help."**



To see all of our body cooling and therapy products, visit [polarproducts.com](http://polarproducts.com).

1.800.763.8423 | [polar@polarproducts.com](mailto:polar@polarproducts.com)

[polarproducts.com](http://polarproducts.com)

